



CUSTOM STABLE FITTING SET MENU

Choose from these clinically proven 24 lens fitting sets (A, B, C or D), or customize your 12 lens set by choosing two 6 Packs out of set A, B, C or D.



14.8 Custom Stable Elite - Oblate

Designed for pediatrics, lower corneal elevations, regular corneas, small aperture/HVID

14.8 Custom Stable Elite - Prolate

Designed for higher cornea elevations, small aperture/HVID

15.8 Custom Stable Elite - Oblate

Designed for regular corneas, high astigmatism, post-surgery (refractive and corneal)

15.8 Custom Stable Elite - Prolate

Designed for all forms of Keratoconus and ectasias

16.8 Custom Stable Elite - Oblate

Designed for larger aperture/HVID, post-surgery (refractive and corneal)

16.8 Custom Stable Elite - Prolate

Designed for larger aperture/HVID, all forms of Keratoconus and ectasias

17.8 Custom Stable Elite - Oblate

Designed for larger aperture/HVID, extreme clearance needs.

17.8 Custom Stable Elite - Prolate

Designed for larger aperture/HVID, extreme clearance needs, all forms of Keratoconus and ectasias.

Set A	6 Pack of	14.8 Custom Stable Elite	Oblate
	6 Pack of	14.8 Custom Stable Elite	Prolate
	6 Pack of	15.8 Custom Stable Elite	Oblate
	6 Pack of	15.8 Custom Stable Elite	Prolate
Set B	6 Pack of	15.8 Custom Stable Elite	Oblate
	6 Pack of	15.8 Custom Stable Elite	Oblate
	6 Pack of	15.8 Custom Stable Elite	Prolate
	6 Pack of	15.8 Custom Stable Elite	Prolate
Set C	6 Pack of	15.8 Custom Stable Elite	Prolate
	6 Pack of	15.8 Custom Stable Elite	Oblate
	6 Pack of	16.8 Custom Stable Elite	Prolate
	6 Pack of	16.8 Custom Stable Elite	Oblate
Set D	6 Pack of	16.8 Custom Stable Elite	Prolate
	6 Pack of	16.8 Custom Stable Elite	Oblate
	6 Pack of	17.8 Custom Stable Elite	Prolate
	6 Pack of	17.8 Custom Stable Elite	Oblate